

Tasohios NW (mm)	Keskihios	Sivukaula (kulma)	Sivukaula (suora)	Kannen korkeus	max. paine (bar)	Tuotenro
<b>NW60</b>						
60	29/32			125	1,5	D4200
60	29/32	1 x 14/23		125	1,5	D4201
60	29/32	2 x 14/23		125	1,5	D4202
60	29/32	1 x 14/23	1 x 14/23	125	1,5	D4203
60	29/32	1 x 19/26		125	1,5	D4204
60	29/32		1 x 19/26	125	1,5	D4205
60	29/32	2 x 14/23	1 x 19/26	125	1,5	D4206
60	29/32	2 x 19/26	1 x 14/23	125	1,5	D4207
<b>NW100</b>						
100	29/32			125	1,0	D4210
100	29/32	1 x 14/23		125	1,0	D4211
100	29/32	2 x 14/23		125	1,0	D4212
100	29/32	3 x 14/23		125	1,0	D4213
100	29/32	1 x 29/32		125	1,0	D4214
100	29/32	2 x 29/32	1 x 14/23	125	1,0	D4215
100	29/32	3 x 29/32		125	1,0	D4216
100	29/32	2 x 29/32	2 x 14/23	125	1,0	D4217
100	45/40			150	1,0	D4218
100	45/40	1 x 29/32		150	1,0	D4219
100	45/40	2 x 29/32	1 x 14/23	150	1,0	D4220
<b>NW120</b>						
120	29/32			125	1,0	D4235
120	29/32	1 x 29/32		125	1,0	D4236
120	29/32	2 x 29/32		125	1,0	D4237
120	29/32	3 x 29/32		125	1,0	D4238
120	29/32	1 x 29/32	1 x 14/23	125	1,0	D4239
120	29/32	2 x 29/32	1 x 14/23	125	1,0	D4240
120	29/32	3 x 29/32	1 x 14/23	125	1,0	D4241
120	29/32	4 x 29/32		125	1,0	D4242
120	29/32	1 x 14/23	1 x 14/23	125	1,0	D4243
120	29/32	2 x 14/23	1 x 14/23	125	1,0	D4244
120	29/32	1 x 29/32	1 x 29/32	125	1,0	D4245
120	29/32	1 x 29/32	1 x 19/26	125	1,0	D4246
120	29/32	2 x 29/32	1 x 19/26	125	1,0	D4247
120	29/32	3 x 29/32	1 x 19/26	125	1,0	D4248
120	45/40	1 x 14/23		150	1,0	D4249
120	45/40	2 x 14/23	1 x 14/23	150	1,0	D4250
120	45/40	2 x 14/23	2 x 29/32	150	1,0	D4251
120	45/40	2 x 29/32		150	1,0	D4252
120	45/40	3 x 29/32	1 x 14/23	150	1,0	D4253
120	45/40	2 x 29/32	2 x 14/23	150	1,0	D4254

<b>NW150</b>						
150	29/32	1 x 29/32		125	1,0	D4261
150	29/32	2 x 29/32		125	1,0	D4262
150	29/32	3 x 29/32		125	1,0	D4263
150	29/32	1 x 29/32	1 x 14/23	125	1,0	D4264
150	29/32	2 x 29/32	1 x 14/23	125	1,0	D4265
150	29/32	3 x 29/32	1 x 14/23	125	1,0	D4266
150	29/32	4 x 29/32		125	1,0	D4267
150	29/32	1 x 14/23	1 x 14/23	125	1,0	D4268
150	29/32	2 x 14/23	1 x 14/23	125	1,0	D4269
150	29/32	1 x 29/32	1 x 29/32	125	1,0	D4270
150	29/32	1 x 29/32	1 x 19/26	125	1,0	D4271
150	29/32	2 x 29/32	1 x 19/26	125	1,0	D4272
150	29/32	3 x 29/32	1 x 19/26	125	1,0	D4273
150	45/40	1 x 14/23		150	1,0	D4274
150	45/40	2 x 14/23	1 x 14/23	150	1,0	D4275
150	45/40	2 x 14/23	2 x 29/32	150	1,0	D4276
150	45/40	2 x 29/32		150	1,0	D4277
150	45/40	3 x 29/32	1 x 14/23	150	1,0	D4278
150	45/40	2 x 29/32	2 x 14/23	150	1,0	D4279
<b>NW200</b>						
200	29/32			145	1,0	D4285
200	29/32	1 x 29/32		145	1,0	D4286
200	29/32	2 x 29/32		145	1,0	D4287
200	29/32	3 x 29/32		145	1,0	D4288
200	29/32	1 x 29/32	1 x 14/23	145	1,0	D4289
200	29/32	2 x 29/32	1 x 14/23	145	1,0	D4290
200	29/32	3 x 29/32	1 x 14/23	145	1,0	D4291
200	29/32	4 x 29/32		145	1,0	D4292
200	29/32	1 x 14/23	1 x 14/23	145	1,0	D4293
200	29/32	2 x 14/23	1 x 14/23	145	1,0	D4294
200	29/32	1 x 29/32	1 x 29/32	145	1,0	D4295
200	29/32	1 x 29/32	1 x 19/26	145	1,0	D4296
200	29/32	2 x 29/32	1 x 19/26	145	1,0	D4297
200	29/32	3 x 29/32	1 x 19/26	145	1,0	D4298
200	45/40	1 x 14/23		170	1,0	D4299
200	45/40	2 x 14/23	1 x 14/23	170	1,0	D4300
200	45/40	2 x 14/23	2 x 29/32	170	1,0	D4301
200	45/40	2 x 29/32		170	1,0	D4302
200	45/40	3 x 29/32	1 x 14/23	170	1,0	D4303
200	45/40	2 x 29/32	2 x 14/23	170	1,0	D4304